

Cochlear™ Baha® Attract System

Cochlear Baha 5 Systems

Cochlear™ Baha® 5 Systems

Introducing the Cochlear Baha 5 Systems



2.4 GHz Wireless



smart App Portfolio



new Fitting Software



new Baha 5 Sound Processor



new Baha Softband



UPGRADE PATH

TRANSITION PATH

new SP Magnet with Colour covers



new 14 mm DermaLock Abutment



stable BI300 Implant



Cochlear™ Baha® 5 Systems

The Cochlear Baha Attract System



2.4 GHz Wireless



smart App Portfolio



new Fitting Software



new Baha 5 Sound Processor



new SP Magnet with
Colour covers



stable BI300 Implant

Cochlear™ Baha® Attract System

An invisible connection



- Discreet and aesthetically appealing
- No need for daily care
- Reduced risk for infections and trauma



Cochlear™ Baha® Attract System

The preferred paediatric choice



When transitioning children from a Baha Softband, the Baha Attract System offers particular benefits:

- Excellent hearing performance
- Good aesthetic outcomes
- Comfort and retention
- Low risk for trauma and infections
- Transition pathway



* In the United States and Canada, the placement of a bone anchored implant is contraindicated in children below the age of five.

Cochlear™ Baha® Attract System

Smart connectivity



The Baha 5 Sound Processor is the only bone conduction sound processor that provides direct-to-device wireless hearing without the need for additional body-worn hardware.



- ✓ Enjoy music
- ✓ Watch movies
- ✓ Stream phone calls
- ✓ Control volume
- ✓ Check device status




Cochlear™ Baha® Attract System

Designed to maximise sound transmission

The Baha Attract System is designed to maximize sound transmission through the skin.

Maximum contact surface area with the **Baha SoftWear™ Pad** for better sound transmission



More powerful and advanced transducer technology



Efficient single point sound transmission

Cochlear™ Baha® Attract System

Unique Baha SoftWear™ Pad




The unique Baha SoftWear™ Pad maximises contact area by adapting to the contours of the skin¹⁻².



Larger contact area leads to:

- Better sound transmission
- Spreading out the magnetic pressure over a larger area enables greater retention force¹⁻²

Largest contact surface area



¹ Flynn M.C. Design concept and technological considerations for the Cochlear Baha 4 Attract System. Cochlear Bone Anchored Solutions AB, Tech brief E82744, 2013.
² Fyrlund, H., Pressure with Attract w/wo SoftWear Pad relative Sophono, Cochlear Bone Anchored Solutions, Technical report no 633480, May 2015


Cochlear™ Baha® Attract System

Better sound transmission



By gently adapting to the contours of the skin, the unique Baha SoftWear Pad maximises the magnet to skin contact area, which improves the sound transmission by an average of **5.4 dB**¹.

Sound transmission with and without Baha SoftWear Pad





- Baha transducer with Baha SoftWear Pad on Artificial Mastoid
- - - Baha transducer without Baha SoftWear Pad on Artificial Mastoid


Cochlear™ Baha® Attract System

Better comfort and retention

Rigid sound processor magnet



Baha Attract System with Baha SoftWear™ Pad




Uneven load distribution¹ with a rigid sound processor plate can cause soft tissue discomfort and complications.

Optimal load distribution¹ with unique Baha SoftWear Pad enables better wearing comfort and less risk of skin irritations.

1.Flynn M.C, Design concept and technological considerations for the Cochlear Baha 4 Attract System. Cochlear Bone Anchored Solutions AB, Tech brief E82744, 2013.

Cochlear™ Baha® Attract System

Efficient single-point sound transmission




MULTIPLE SCREWS



SINGLE-POINT

Multiple screw systems spread vibration energy over several contact points.

By focusing vibrations in one single point, the Baha Attract System is designed to transmit sound more efficiently to the inner ear.¹




Most stable Baha implant ever.²

1. Flynn M.C. Design concept and technological considerations for the Cochlear Baha 4 Attract System. Cochlear Bone Anchored Solutions AB, Tech brief E82744, 2013.

2. Nelissen RC, Stalfors J, de Wolf MJ, Flynn MC, Wigren S, Eeg-Olofsson M, Green K, Rothera MP, Mylanus EA, & Hol MK. Long-term stability, survival, and tolerability of a novel osseointegrated implant for bone conduction hearing: 3-year data from a multicenter, randomized, controlled, clinical investigation. Otol Neurotol. 2014, 35(8): 1486-91.

Cochlear™ Baha® Attract System

Choice of magnets for individual patients needs



New rounded edges
for maximum discretion

Same comprehensive range of
strengths to suit individual
patients


SP Magnet	Soft tissue thickness, mm			
	3	4	5	6
No. 6				Extra
No. 5			Extra	Normal
No. 4		Extra	Normal	Start
No. 3	Extra	Normal	Start	Low
No. 2	Normal	Start	Low	Low
No. 1	Start	Low	Low	Low
No. 2 with 2 SoftWear Pads	Low			
No. 1 with 2 SoftWear Pads	Low			

Cochlear™ Baha® Attract System


Smooth magnet attached to a strong foundation



Sleek titanium casing with rounded edges to minimise soft tissue irritation.



Established coupling to facilitate removal of the implant magnet for future upgrades and possible MRI > 1.5 Tesla.



Strong foundation for efficient single-point sound transmission and future hearing options.

Cochlear™ Baha® Attract System

Straightforward, proven and reliable



- **Proven procedure**
 - More than 5000 patients implanted.
 - 0% implant loss reported to Cochlear (June 2015).**
- Short surgical time:
 - Surgery times down to 24 minutes have been published¹.
- Uneventful post-operative healing enables fitting in four weeks¹.




***Figure based on reported implant loss compared with total sales*

1. Briggs R, Van Hasselt A, Luntz M, Goycoolea M, Wigren S, Weber P, Smeds H, Flynn M, Cowan R. Clinical performance of a new magnetic bone conduction hearing implant system: results from a prospective, multicenter, clinical investigation. *Otol Neurotol.* 2015;36(5):834-41.

Less invasive surgery and options for life

1. Less invasive surgery

The implant and magnet can be placed without the need for extensive bone removal, dura exposure or preoperative CT scanning.



2. Efficient transmission

The osseointegrated implant provides efficient single point sound transmission.

3. Upgrade options

Patients can upgrade their hearing to a more powerful sound processor or by transition to direct bone conduction.



Baha Attract System Summary



SP Magnets and Colour covers designed to blend in

Smallest sound processor in the industry

Unique Baha SoftWear™ Pad for superior comfort and better sound transmission


A portfolio of sound processors

Transition path to an abutment

True wireless accessories and Made for Iphone connectivity with most complete App portfolio.

System components designed for most efficient sound transmission

Cochlear support with life-time commitment





Cochlear®

Hear now. And always